




CompuMed Exam Website

Uploading via Web-Browser (no software install required)


1. To upload directly from the website, go to the **[My Studies]** tab at the top of the screen, then find the **[Study Actions]** tab toward the bottom of the screen below the exam list.




2. Click on the **[Study Actions]** tab, scroll up to find the **[Study Upload]** button.




3. In the **[Upload Files to Server]** window, click the **[Choose Files]** options to browse files.




4. Navigate to the drive with the **DICOM files** (usually the CD/DVD drive).



5. Select all the files in the folder, by clicking and dragging or ctrl + A, and click **[Open]**.



Most DICOM disks will have multiple folders.



Search through the folders until you find the DICOM files. Here's an example of what DICOM files generally look like on a disc.

6. On the **[Upload Files to Server]** window, you'll see the number of files select
7. Click on **[Upload Files]** to start the upload.

If you have any questions, please call 310-258-5012 or support@compumedinc.com

www.compumedinc.com